




Recombinant Junin mammarenavirus Nucleoprotein (N)

Product Code	CSB-EP318492JAM
Relevance	Encapsidates the genome, protecting it from nucleases. The encapsidated genomic RNA is termed the nucleocapsid (NC). Serves as template for viral transcription and replication. The increased presence of protein N in host cell does not seem to trigger the switch from transcription to replication as observed in other negative strain RNA viruses. Disables the host innate defense by interfering with beta interferon (IFNB) production through the inhibition of host IRF3 phosphorylation and activation by host IKBKE. Through the interaction with host IKBKE, strongly inhibits the phosphorylation and nuclear translocation of IRF3, a protein involved in IFN activation pathway, leading to the inhibition of IFNB and IRF3-dependent promoters activation (By similarity). Interacts with host IKBKE (via Protein kinase domain); the interaction inhibits IKBKE kinase activity.
Abbreviation	Recombinant Junin mammarenavirus N protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P14239
Alias	Nucleocapsid protein Protein N
Product Type	Recombinant Protein
Immunogen Species	Junin mammarenavirus (JUNV) (Junn mammarenavirus)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MAHSKEVPSFRWTQSLRRGLSQFTQTVKSDVLKDAKLIADSIDFNQVAQVQRA LRKTKRGEEDLNKLRDLNKEVDRLMSMRSVQRNTVFKAGDLGRVERMELAS GLGNLKTFRRAETGSQGVYMGNLSQSQLAKRSEILRTLGFQQQGTGGNGV VRVWDVKDPSKLNQFGSVPALTIACMTVQGGETMNSVIQALTSLGLLYTVKY PNLSLDLRLTQEHDCQLQIVTKDESSINISGYNFSLSAAVKAGASILDDGNMLETI RVTPDNFSSLIKSTIQVKRREGMFIDEKPGNRNPYENLLYKLCLSGDGWVPYIGS RSQIIGRSWDNTSIDLTRKPVAGPRQPEKNGQNLRLANLTEIQEAVIREAVGKL DPTNTLWLDIEGPATDPVEMALFQPAGSKYIHCFRKPHDEKGFKNRSRSHGI LMKDIEDAMPGVLSYVIGLLPPDMVVTTQGSDDIRKLFDLHGRRDLKLVDVRLT SEQARQFDQQVWEKFGHLCKHHNGVVVSKKKRDKDAPFKLASSEPHCALLD CIMFQSVLDGKLYEEELTPLLPPSLLFLPKAAYAL
Research Area	Others
Source	E.coli
Target Names	N
Expression Region	1-564aa



Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	79.0kDa
Protein Length	Full Length

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-----------------------	---

Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
-------------------	---