





Recombinant Streptomyces alboniger Puromycin N-acetyltransferase (pac)

Product Code	CSB-EP319512SNI
Relevance	Detoxification of puromycin.
Abbreviation	Recombinant Streptomyces alboniger pac protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P13249
Product Type	Recombinant Protein
Immunogen Species	Streptomyces alboniger
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MTEYKPTVRLATRDDVPRAVRTLAAAFADYPATRHTVDPDRHIERVTELQELF LTRVGLDIGKVWVADDGAAVAVWTTTPESVEAGAVFAEIGPRMAELSGSRLAA QQQMEGLLAPHRPKPAWFLATVGVSPDHQKGKGLGSAVVLPGVEAAERAGV PAFLETSAPRNLPFYERLGFTVTADVEVPEGPRTWCMTRKPGA
Source	E.coli
Target Names	pac
Expression Region	1-199aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	37.5kDa
Protein Length	Full Length


Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP319512SNI could indicate that this peptide derived from E.coli-expressed *Streptomyces alboniger* pac.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP319512SNI could indicate that this peptide derived from E.coli-expressed *Streptomyces alboniger* pac.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.