




# Recombinant Moraxella bovis Fimbrial protein 1

<b>Product Code</b>	CSB-EP322675MHZ
<b>Abbreviation</b>	Recombinant Moraxella bovis Fimbrial protein 1 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P20657
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Moraxella bovis
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	FTLIELMIVIAIIGILAAIALPAYQDYISKSQTTRVSGELAAGKTAVDAAALFEGKTPV LSEESSTSKENIGLTSSETSTKPRSNLMAVELTGFADNGAGTISATLGNKANK DIAKTIVITQERTTDGVWTCKIDGSQAQYKEKFNPTGCVKK
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	N/A
<b>Protein Names</b>	Recommended name: Fimbrial protein 1Alternative name(s): Alpha-pilin Fimbrial protein I I pilin
<b>Expression Region</b>	7-159aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	29.0 kDa
<b>Protein Length</b>	Full Length of Mature Protein

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



## **Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.