




Recombinant Escherichia coli Thiosulfate-binding protein (cysP)

Product Code	CSB-EP325847ENV
Abbreviation	Recombinant E.coli cysP protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P16700
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Escherichia coli (strain K12)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	TELLNSSYDVSRELFAALNPPFEQQWAKDNGGDKLTIKQSHAGSSKQALAILQ GLKADVVTYNQVTDVQILHDKGKLI PADWQSRLPNNSSPFYSTMGFLVRKGNP KNIHDWNDLVRSDVKLIFPNPKTSGNARYTYLAAWGAADKADGGDKGKTEQF MTQFLKNVEVFD TGGRGATTTFAERGLGDVLISFESEVNNIRKQYEAQGFEVVI PKTNILAEFPVAWVDKNVQANGTEKAAKAYLNWLYSPQAQTIITDYYYYRVNNP EVM DKLKDKFPQTELF RVEDKFGSWPEVMKTHFTSGGELDKLLAAGR N
Research Area	Others
Source	E.coli
Target Names	cysP
Protein Names	Recommended name: Thiosulfate-binding protein
Expression Region	26-338aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	48.0 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.