



# Recombinant Severe acute respiratory syndrome coronavirus 2 Nucleoprotein (N) (P13L,Δ31-33,E136D,R203K,G204R,S413R), partial

<b>Product Code</b>	CSB-EP3325GMY(M19)
<b>Abbreviation</b>	Recombinant SARS-CoV-2 N protein (P13L,Δ31-33,E136D,R203K,G204R,S413R), partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0DTC9
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
<b>Product Type</b>	Recombinant Severe acute respiratory syndrome coronavirus 2 Nucleoprotein(N)(P13L,Δ31-33,E136D,R203K,G204R,S413R),partial
<b>Immunogen Species</b>	Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	MSDNGPQNQRNAPRITFGGSDSTGNSQNGGARSKQRRPQGLPNNTASWF TALTQHGKEDLKFPRGQGVPIINTNSSPDDQIGYYRRATRRIRGGDGKMKDLS PRWYFYLLGTGPEAGLPYGANKDGIWVATDGALNTPKDHIGTRNPANNAIV LQLPQGTTLPKGFYAEGSRGGSQASSRSSSRNSSRNSTPGSSKRTSPAR MAGNGGDAALALLLDRLNQLESKMSGKGGQQGQTVTKKSAAEASKKPRQ KRTATKAYNVTQAFGRRGPEQTQGNFGDQELIRQGTDYKHWPQIAQFAPSAS AFFGMSRIGMEVTPSGTWLTYTGAIKLDDKDPNFKDQVILLNKHIDAYKTFPPT EPKDKKKKKADETQALPQRQKKQQTVLLPAADLDDFSKQLQQSMSRADSTQ A
<b>Research Area</b>	Microbiology
<b>Source</b>	E.coli
<b>Target Names</b>	N
<b>Expression Region</b>	1-419aa(P13L,Δ31-33,E136D,R203K,G204R,S413R)
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged



---


**Mol. Weight** 49.5 kDa

---

**Protein Length** Partial

---

**Image**



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.