




Recombinant *Saccharomyces cerevisiae* A-agglutinin-binding subunit (AGA2)

Product Code	CSB-EP335766SVG
Relevance	Receptor binding subunit of the a-agglutinin heterodimer. <i>S.cerevisiae</i> a and alpha cells express the complementary cell surface glycoproteins a-agglutinin and alpha-agglutinin, respectively, which interact with one another to promote cellular aggregation during mating.
Abbreviation	Recombinant <i>Saccharomyces cerevisiae</i> AGA2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32781
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	QELTTICEQIPSPITLESTPYSLSTTTILANGKAMQGVFEYYKSVTFVSNCGSHP STTSKGSPINTQYVF
Research Area	Others
Source	<i>E.coli</i>
Target Names	AGA2
Protein Names	Recommended name: A-agglutinin-binding subunit
Expression Region	19-87aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	23.5kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.