




Recombinant Human Trifunctional enzyme subunit alpha, mitochondrial (HADHA)

Product Code	CSB-EP340324HU
Abbreviation	Recombinant Human HADHA protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40939
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Trifunctional enzyme subunit alpha, mitochondrial (HADHA)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	TRTHINYGVKGDVAVVRINSPNSKVNTLSKELHSEFSEVMNEIWASDQIRSAVL ISSKPGCFIAGADINMLAACKTLQEVTLQSLQEAQRIVEKLEKSTKPIVAAINGSCL GGGLEVAISCQYRIATKDRKTVLGTPEVLLGALPGAGGTQRLPKMVGVPAAALD MMLTGRSIRADRAKKMGLVDQLVEPLGPGLKPPEERTIEYLEEVAITFAKGLAD KKISPKRDKGLVEKLTAYAMTIPFVRQQVYKKVEEKVRKQTKGLYPAPLKIIDVV KTGIEQGS DAGYLCE SQKFGE LVMTKESKALMGLYHGQVLCKKNKFGAPQKD VKHLAILGAGLMGAGIAQVSVDKGLKTLKDATLTALDRGQQQVFKGLNDKVKK KALTSEFERDSIFSNTLGQLDYQGF EKADMVIEAVFEDLSLKHRVLKEVEAVIPD HCIFASNTSALPISEIAAVSKRPEKVIGMHYFSPVDKMQLLEITTEKTSKDTSAS AVAVGLKQKGKVIIVVKDGP GFYTTRCLAPMMSEVIRILQEGVDPKKLDSLTTSF GFPVGAATLVDEVGVDVAKHVAEDLGKVFGERFGGGNPELLTQMVSKGFLGR KSGKGFYIYQEGVKRKDLNSDMSILASLKLPKSEVSSDEDIQFRLVTRFVNE AVMCLQEGILATPAEGDIGAVFGLGFPPCLGGPFRFVDLYGAQKIVDRLKKYEA AYGKQFTPCQLLADHANS PNKKFYQ
Research Area	Cardiovascular
Source	E.coli
Target Names	HADHA
Expression Region	37-763aa
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	86.0 kDa


Protein Length

Full Length of Mature Protein

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.