




# Recombinant Streptococcus mutans serotype c Cell surface antigen I/II (spaP), partial

<b>Product Code</b>	CSB-EP340548FMR
<b>Abbreviation</b>	Recombinant Streptococcus mutans serotype c spaP protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P23504
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
<b>Product Type</b>	Recombinant Streptococcus mutans serotype c Cell surface antigen I/II(spap),partial
<b>Immunogen Species</b>	Streptococcus mutans serotype c (strain ATCC 700610 / UA159)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	PQVNKEIRNNNDINIDRTLVAKQSVVKFQLKTADLPAGRDETTSFVLVDPLPSG YQFNPEATKAASPGFDVTDNATNTVTFKATAATLATFNADLTKSVATIYPTVV GQVLNDGATYKNNFTLTVNDAYGIKSNVVRVTTGPKNDPDPNPNNNYIKPTKV NKNENGVIDGKTVLGASTNYYELTWDLDQYKNDRSSADTIQKGFYYVDDYP EEALELRQDLVKITDANGNEVTGVSVDNYTNLEAAPQEIRDVLSKAGIRPKGAF QIFRADNPREFYDITYVKTGIDLKIVSPMVVKKQMGGTGGSYENQAYQIDFGNG YASNIINNVPKINPKKDVTLLDPADTNNDGQTIPLNTVFNYRLIGGIIPADHSE ELFEYNFYDDYDQTDGHTGQYKVFQVDITFKDGSIIKSGAELTQYTAEVDT AKGAIKFKFAFLRSVSIDSAFQAESYIQMKRIAVGTFENTYINTVNGVTYSSNT VKTTTPEDPTDPTDPQDPSSPRTSTVINYPKQSTAYQPSSVQETLPNT
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	spaP
<b>Expression Region</b>	998-1532aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	C-terminal 6xHis-tagged
<b>Mol. Weight</b>	60.0 kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .