



Recombinant *Eschscholzia californica* Reticuline oxidase (BBE1)

Product Code	CSB-EP341191EOT
Relevance	Essential to the formation of benzophenanthridine alkaloids in the response of plants to pathogenic attack. Catalyzes the stereospecific conversion of the N-methyl moiety of (S)-reticuline into the berberine bridge carbon of (S)-scoulerine.
Abbreviation	Recombinant <i>Eschscholzia californica</i> BBE1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P30986
Alias	Berberine bridge-forming enzyme ;BBETetrahydroprotoberberine synthase
Product Type	Recombinant Protein
Immunogen Species	<i>Eschscholzia californica</i> (California poppy)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	GNDLLSCLTFNGVRNHTVFSADSDSDFNRFLHLSIQNPLFQNSLISKPSAILPG SKEELSENTIRCIRKGSWTIRLRSGGHSYEGLSYTS DTPFILIDLMNLNRV SIDLES ETAWVESGSTLGELYAITESSSKLGFTAGWCPTVGTGGHISGGGFGMM SRK YGLAADNVVDAILIDANGAILDRQAMGEDVFWAIRGGGGGVWGAIYAWKIKLL PVPEKVTVFRVTKNVAIDEATSL LHKWQFVAEELEEDFTLSVLGGADEKQVWL TMLGFHFGLKTVAKSTFDLLFPELGLVEEDYLEMSWGESFAYLAGLETVSQ L N NRFLKFDERAFKTKVDLTKEPLPSKAFYGLLERLSKEPNGFIALNGFGGQMSKI SSDFTPFPHRSGTRLMVEYIVAWNQSEQKKKTEFLDWLEKVYEFMKPFVSKN PRLGYVNHIDLDLGGIDWGNKTVVNNAIEISR SWGESYFLSNYERLIRAKTLIDP NNVFNHPQSIPPMANFDYLEKTLGSDGGEVVI
Research Area	Others
Source	E.coli
Target Names	BBE1
Expression Region	24-538aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	61.4kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.