




Recombinant Leiurus hebraeus Alpha-like toxin Lqh7

Product Code	CSB-EP348581LDS
Abbreviation	Recombinant Leiurus hebraeus Alpha-like toxin Lqh7 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P59357
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Leiurus hebraeus (Deathstalker scorpion) (Leiurus quinquestriatus hebraeus)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	VRDGYIAKPENCAHHCFPGSSGCDTLCKENGGTGGHCGFKVGHGTACWCNALPDKVGIIVDGVKCH
Research Area	Others
Source	E.coli
Target Names	N/A
Protein Names	Recommended name: Alpha-like toxin Lqh7Alternative name(s): Lqh VII Short name= LqhVII
Expression Region	1-66aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	11.8 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.