




Recombinant Brucella melitensis biotype 1 26 kDa periplasmic immunogenic protein (bp26)

Product Code	CSB-EP358424BMO
Abbreviation	Recombinant Brucella melitensis biotype 1 bp26 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A3U8
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Brucella melitensis biotype 1 (strain 16M / ATCC 23456 / NCTC 10094)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	QENQMTTQPARIAVTGEGMMTASPDMAILNLSVLRQAKTAREAMTANNEAMT KVLDMKKAGIEDRDLQTGGINIQPIYVYPDDKNNLKEPTITGYSVSTSLTVRVR ELANVGKILDESVTLGVNQGGDLNLVNDNPSAVINEARKRAVANAIKAKTLAD AAGVGLGRVVEISELSRPPMPMPIARGQFRTMLAAAPDNSVPIAAGENSYNVS VNVVFEIK
Research Area	Others
Source	E.coli
Target Names	bp26
Expression Region	29-250aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	31.1 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.