




Recombinant Human papillomavirus type 11

Regulatory protein E2 (E2)

Product Code	CSB-EP361192HMG
Abbreviation	Recombinant Human papillomavirus 11 E2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04015
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Human papillomavirus type 11
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MEAIKRLDACQDQLLELYEENSIDIHKHIMHWKCI RLESVLLHKAKQMGLSHIG LQVVPPLTVSETKGHNAIEMQMHLESLAKTQYGV EPWTLQDTSYEMWLTTPPK RCFKKQGNTVEVKFDGCEDNVMEYVWVTHIYLQDNDSWVKVTSSVDAKGIYY TCGQFKTYVNFNKEAQKYGSTNHWEVCYGSTVICSPASVSSTVREVSIAEPT TYTPAQTTAPTVSACTTEDGVSAPPRKRARGPSTNNTLCVANIRSVDSTINNIV TDNYNKHQRRNNCHSAATPIVQLQGDSNCLKCFRYRLNDKYKHLFELASSTW HWASPEAPHKNAIVTLTYSSEEQRQQFLNSVKIPPTIRHKVGFMSLHLL
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	E2
Expression Region	1-367aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	49.2 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Description

The recombinant Human papillomavirus type 11 Regulatory protein E2 expression involves transfecting e.coli cells with a DNA expression vector containing the gene encoding the recombinant protein of interest (1-367aa) and then culturing the cells so that they express the desired protein. The protein features a N-terminal 10xHis tag and C-terminal Myc tag fusion. The recombinant Human papillomavirus type 11 Regulatory protein E2 is collected and purified from the cell lysate through affinity purification. Its purity is greater than 85%, as measured by SDS-PAGE.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.