




Recombinant Yellow fever virus Genome polyprotein, partial, Biotinylated

Product Code	CSB-EP365905YAC1g10-B
Abbreviation	Recombinant Yellow fever virus Genome polyprotein, partial, Biotinylated
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03314
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Yellow fever virus Genome polyprotein, partial, Biotinylated
Immunogen Species	Yellow fever virus (strain 17D vaccine) (YFV)
Sensitivity	Not Test
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	AHCIGITDRDFIEGVHGGTWVSATLEQDKCVTVMAPDKPSLDISLETVAIDRPA EVRKVCYNAVLTHVKINDKCPSTGEAHLAEENEGDNACKRTYSDRGWGNGC GLFGKGSIVACAKFTCAKSMSLFEVDQTKIQYVIRAQLHVGAKQENWNTDIKTL KFDALSGSQEVEFIGYGKATLECQVQTAVDGNSYIAEMETESWIVDRQWAQ DLTLPWQSGSGGVWREMHHLVEFEPHAATIRVLALGNQEGSLKTALTGAMR VTKDTNDNNLYKLHGGHVSCRVKLSALTCLKGTSYKICTDKMFFVKNPTDTGHG TVVMQVKVSKGAPCRIPVIVADDLTAANKGILVTVNPIASTNDDEVILIEVNPPFG DSYIIVGRGDSRLTYQWHKEGSSIGKLFTQTMKGVERLAVMGDTAWDFSSAG GFFTSVGKGIHTVFGSAFQGL
Research Area	Others
Source	E.coli
Target Names	N/A
Expression Region	286-730aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal MBP-tagged and C-terminal 6xHis-Avi-tagged
Mol. Weight	96.2 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.