




Recombinant Rabies virus Glycoprotein (G), partial

Product Code	CSB-EP365970RIDa0
Abbreviation	Recombinant Rabies virus Glycoprotein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03524
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Rabies virus (strain ERA)(RABV)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	KFPIYTILDKLGPWSPIDIHHLSCPNNLVVEDEGCTNLSGFSYMELEKVGILAIK MNGFTCTGVVTEAETYTNFVGIVTTTFKRKHFRPTPDACRAAYNWKMAGDP RYEESLHNPYPDYRWLRTVKTTKESLVIISPSVADLDPYDRSLHSRVFSPGKCS GVAVSSTYCSTNHDIWMPENPRLGMSCDIFTNSRGKRASKGSETCGFVDE RGLYKSLKGACKLKLGGVLGLRLMDGTWVAMQTSNETKWCPPDQLVNLHDF RSDEIEHLVVEELVRKREECLDALESIMTTKSVSFRRLSHLRKLVPFGFKAYTIF NKTLMEDIAHYKSVRTWNEILPSKGCLRGGGRCHPHVNGVFFNGIILGPDGNV LIPEMQSSLLQQHMELLESSVIPLVHPLADPSTVFKDGDGAEDFVEVHLPDVHN QVSGVDLGLPNWGKY
Research Area	Signal Transduction
Source	E.coli
Target Names	G
Expression Region	20-459aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	55.0 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.