



Recombinant Plasmodium ovale Tryptophan-rich protein (PowCR01_000150100), partial

Product Code	CSB-EP3733GPG
Abbreviation	Recombinant Plasmodium ovale PowCR01_000150100 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0A1C3KIR3
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Plasmodium ovale (malaria parasite P. ovale)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	ASNEKNLLLLDVENDPLSFLKPITDTANPINTANPINMANPINMAKPINTITPINAI TPINAITPINAITPIKTITPINTITSINTITPMNTITPISTTKPVNTTNTINTTNTLATLD TIAKPTTQAEIKQHMGAAWAKQIHDEWEIFVETRRQEVENSEKEKEQEWDN WLTELQPRWENFRSFVDKKKEKWKIEKKELEWGKWIKAMEGKWMTYKENIDK RYFPGNMKDAAADSEWVNNIKEEVQKKMIIDLKKWISGNNINLYKWILRDWND WKKKRLDEWNNREWKCEEDKYWENYSRKGNGWDPLNILLKKKCEKWIERKK REKDHWN DLTRRL EITYLSKKHMEWEKWKDDKCEWYSKWIMFF
Research Area	Neuroscience
Source	E.coli
Target Names	PowCR01_000150100
Expression Region	32-399aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	51.5 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.