



Recombinant Human Striatin-3 (STRN3)


Product Code	CSB-EP614388HUd7
Abbreviation	Recombinant Human STRN3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
Uniprot No.	Q13033
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Human Striatin-3 (STRN3)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	<p> MDELAGGGGGPGMAAPPRQQQGPGGNLGLSPGGNGAAGGGGPPASEGA GPAAGPELSRPQQYTIPGILHYIQHEWARFEMERAHWEVERAELQARIAFLQG ERKGQENLKKDLVRRIMLEYALKQERAKYHKLKYGTELNQGDLMPTFESEE TKDTEAPTAPQNSQLTWKQGRQLLRQYLQEVGYTDITLDVRSQRVRSLLGLS NSEPNGSVETKNLEQILNGGESPKQKQGEIKRSSGDVLETNFLENADDSDED EENDMIEGIPEGKDKHRMKNKHKIGNEGLAADLTDDPDTEEALKEFDLFTAED GEGAGEARSSGDGTEWDKDDLSPTAEVWDVDQGLISKLKEQYKKERKGKKG VKRANRTKLYDMIADLGDELPHIPSGIINQSRSASTRMTDHEGARAEAEAPIT FPSGGGKSFIGSDDVLLSVLGLGLADLTVTNDADYSYDLPANKDAFRKTW NPKYTLRSHFDGVRALAFHPVEPVLTASEDHTLKLWNLQKTVPAKKSASLDV EPIYTFRAHIGPVLSLAISSNGEQCFSGGIDATIQQWWNMPSPSVDPTYDYEPNV LAGTLVGHTDAVWGLAYSGIKNQLLSCSADGTVRLWNPQEKLPCTYNGDK KHGIPTSVDFIGCDPAHMTSFNTGSAVIYDLETSQSLVILSSQVDSGLQSNHHI NRVVSHTPLPVTITAHEDRHIKFFDNKTGKMIHSMVAHLDAVTS LAVDPNGIYL MSGSHDCSIRLWNLD SKTCVQEITAHRRKKLDESIYD VAFHSSKAYIASAGADAL AKVFV </p>
Research Area	Others
Source	E.coli
Target Names	STRN3
Expression Region	1-797aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	C-terminal 10xHistagged



Mol. Weight 97.0 kDa

Protein Length Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.