




Recombinant Human Embryonic stem cell-related gene protein (HESRG)

Product Code	CSB-EP632275HU
Abbreviation	Recombinant Human ESRG protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1W209
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MTLFSDSARLHPGEINSLVAHTKPVWWSLHTDAHEIWCRDSDRGTSLGSRIPC PPALCSVRKIHLRPQVLRPTSPRNISPISNPVSGLFLLCSPTSLTIPQPLSPFNLG ATLQSLPSLNFNSFHSLVETKETCFIREPKTPAPVTDWEGSLPLVFNHCRDASL ISRFRPRRDACLGPSPLAASPAFLGQGQVPLNPFSFTLSGKSRFSGAGASTPQ PLLLHP
Research Area	Others
Source	E.coli
Target Names	ESRG
Protein Names	Recommended name: Embryonic stem cell-related gene protein Short name= hES cell-related gene protein
Expression Region	1-222aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	28.2kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.