




Recombinant Human Inactive N-acetylated-alpha-linked acidic dipeptidase-like protein 2 (NAALADL2), partial

Product Code	CSB-EP704507HUa0
Abbreviation	Recombinant Human NAALADL2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q58DX5
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	GYYVHTNCPSDAPSSGTVDPQLYQEILKTIQAEDIKKSFRNLVQLYKNEDDMEI SKKIKTQWTSLGLEDVQFVNYSVLLDLPGPSPSTVTLSSSGQCFHPNGQPCSE EARKDSSQDLLYSYAAYSAGGTLKAEVIDVSYGMADDLKRIRKIKNVTNQIALLK LGKLPLLYKLSSLEKAGFGGVLLYIDPCDLPKTVNPSHDTFMVSLNPGGDPSTP GYPSVDESFRQSRNLTSLLVQPISAPLVAKLISSPKARTKNEACSSLELPNNEI RVVSMQVQTVTKLKTVTNVVGFVMGLTSPDRYIIVGSHHHTAHSYNGQEWAS STAITAFIRALMSKVKRGWRPDRTIVFCSWGGAFTGNIGSYEWGEDFKKVLQ KNVVAYISLHSPIRGNSSLYPVASPSLQQLVVEKNNFNCTRRACPETNISSIQI QGDADYFINHLGVPIVQFAYEDIKTLEGPSFLSEARFSTRATKIEEMDPSFNLHE TITKLSGEVILQIANEPVLPFNALDIALEVQNNLKGDPNTHQLLAMALRLRESA ELFQSDEMRPANDPKERAPIRIRMLNDILQDMEKSFLVKQAPPGFYRNILYHLD EKTSRFSILIEAWEHCKPLASNETLQEALSEVLNSINSAQVYFKAGLDVFKSVLD GKN
Research Area	Cell Biology
Source	E.coli
Target Names	NAALADL2
Expression Region	143-795aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	78.6 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.