




Recombinant Rat 17-beta-hydroxysteroid dehydrogenase 13 (Hsd17b13)

Product Code	CSB-EP711723RA
Abbreviation	Recombinant Rat Hsd17b13 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5M875
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Rattus norvegicus (Rat)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	ESLVKFFIPQRRKSVAGQTVLITGAGHGIGRLTAYEFAKQKSRLVLWDISKHGV EETAACKCRKLGAVVHV FVVD C SNRAEIYKSVDQVKKEVGDI EIVVNNAGAIYPA DLLSTKDEEITKT FEVNILGHFWI IKALLPSMLRRNSGHIVTVASVCGHRVIPYLIP YCSSKFAAVG FHRALTAELDTLGKTGIKTSCLCPVFVNTGFTKNPSTRLWPVLE PDEVARSLIDGILT NKKMIFVPSYINISLIVEMFFPERVLKAINRIQNIQFEAIVGHR TKRK
Research Area	Others
Source	E.coli
Target Names	Hsd17b13
Protein Names	Recommended name: 17-beta-hydroxysteroid dehydrogenase 13 Short name= 17-beta-HSD 13 EC= 1.1.-.- Alternative name(s): Short-chain dehydrogenase/reductase 9
Expression Region	20-300aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	37.4 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.