




Recombinant Mouse Killer cell lectin-like receptor 3 (Klra3), partial

Product Code	CSB-EP720664MO
Relevance	Receptor on natural killer (NK) cells for class I MHC.
Abbreviation	Recombinant Mouse Klra3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64329
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	QYNQHKQEINETLNHHHNCNMQRAFNLKEEMLTNKSIDCRPSNETLEYIKRE QDRWDSKTKTVLDSSRDTGRGVKYWFCYSTKCYFIMNKTTWSGCKANCQH YSVPILKIEDEDELFQRLQHVIPENYWIGLSYDKKKKEWAWIDNGPSKLDMKIR KMNFKSRGCVFLSKARIEDIDCNIPYYCICGKKLDKFPD
Research Area	Others
Source	E.coli
Target Names	Klra3
Protein Names	5E6 Lymphocyte antigen 49c Short name: Ly-49c Nk2.1 T-cell surface glycoprotein Ly-49C Ly-49c, Ly49C
Expression Region	70-266aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	27.6 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.