




Recombinant Mouse E3 ubiquitin-protein ligase TRIM21 (Trim21)

Product Code	CSB-EP733816MO
Relevance	E3 ubiquitin-protein ligase whose activity is dependent on E2 enzymes, UBE2D1, UBE2D2, UBE2E1 and UBE2E2. Forms a ubiquitin ligase complex in cooperation with the E2 UBE2D2 that is used not only for the ubiquitination of USP4 and IKBKB but also for its self-ubiquitination. Component of cullin-RING-based SCF (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complexes such as SCF(SKP2)-like complexes. A TRIM21-containing SCF(SKP2)-like complex is shown to mediate ubiquitination of CDKN1B ('Thr-187' phosphorylated-form), thereby promoting its degradation by the proteasome. Monoubiquitinates IKBKB that will negatively regulates Tax-induced NF-kappa-B signaling. Negatively regulates IFN-beta production post-pathogen recognition by polyubiquitin-mediated degradation of IRF3. Mediates the ubiquitin-mediated proteasomal degradation of IgG1 heavy chain, which is linked to the VCP-mediated ER-associated degradation (ERAD) pathway. Promotes IRF8 ubiquitination, which enhanced the ability of IRF8 to stimulate cytokine genes transcription in macrophages. Plays a role in the regulation of the cell cycle progression. Enhances the decapping activity of DCP2. Exists as a ribonucleoprotein particle present in all mammalian cells studied and composed of a single polypeptide and one of four small RNA molecules. At least two isoforms are present in nucleated and red blood cells, and tissue specific differences in RO/SSA proteins have been identified. The common feature of these proteins is their ability to bind HY RNAs. ²
Abbreviation	Recombinant Mouse Trim21 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62191
Alias	52 kDa Ro protein; 52 kDa ribonucleoprotein autoantigen Ro/SS-A; Ro(SS-A); Sjogren syndrome type A antigen ; SS-A; Tripartite motif-containing protein 21
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MSPSTTSKMSLEKMWEEVTCISICLDPMVEPMSIECGHCFCKECEFVGKNGG SSCPECRQQFLLRNLRPNRHIANMVENLKQIAQNTTKKSTQETHCMKHGKELHL FCEEDGQALCWVCAQSGKHRDHTRVPIEEAAKVYQEKIHVVLEKLRKGKELAE KMEMDLTMQRDWDKRNIDTQKSRIHAEFALQNSLLAQEEQRQLQRLEKDQRE YLRLLGKKEAEAEKNQALQELISELERRIRGSELELLQEVRIILERSGSWNLD LDIDAPDLTSTCPVPGRRKMLRTCWVHITLDRNTANSWLIISKDRRQVRMGDT HQNVSDNKERFSNYPMVLGAQRFSSGKMYWEVDVTQKEAWDLGVCRDSVQ



RKGQFSLSPENGFWTIWLWQDSYEAGTSPQTTLHIQVPPCQIGIFVDYEAGVV
SFYNITDHGSLIYTFSECVFAGPLRPFFNVGFNYSGGNAAPLKLCPLKM

Research Area	Others
Source	E.coli
Target Names	Trim21
Expression Region	1-470aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	70.2kDa
Protein Length	Full Length

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.