




# Recombinant Human Dual specificity protein phosphatase 13 isoform MDSP (DUSP13)

<b>Product Code</b>	CSB-EP734923HU
<b>Relevance</b>	Probable protein tyrosine phosphatase. Has a phosphatase activity-independent regulatory role in MAP3K5/ASK1-mediated apoptosis, preventing MAP3K5/ASK1 inhibition by AKT1. Shows no phosphatase activity on MAPK1/ERK2, MAPK8/JNK, MAPK14/p38 and MAP3K5/ASK1.
<b>Abbreviation</b>	Recombinant Human DUSP13 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6B8I1
<b>Alias</b>	Branching-enzyme interacting DSPMuscle-restricted DSP ;MDSP
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	MAETSLPELGGEDKATPCPSILELEELLRAGKSSCSRVDVWPNLFIGDAATA NNRFELWKLGITHVLNAAHKGLYCQGGPDFYGSSVSYLGVPALDLPDFDISAY FSSAADFIHRA LNTPGAKVLVHCVVGVSR SATLVLAYLMLHQRLSLRQAVITVR QHRWVFPN RGF LHLQ LCR LDQQLRGAGQS
<b>Research Area</b>	Signal Transduction
<b>Source</b>	E.coli
<b>Target Names</b>	DUSP13
<b>Expression Region</b>	1-198aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	36.7kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.