




Recombinant Human BH3-like motif-containing cell death inducer (BLID)

Product Code	CSB-EP818276HU
Relevance	Functions as a proapoptotic molecule through the caspase-dependent mitochondrial pathway of cell death.
Abbreviation	Recombinant Human BLID protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IZY5
Alias	Breast cancer cell protein 2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MVTLLPIEGQEIHFFEILESECVLYTGWIERASGSSIYPEAKARLPLEALLGSNKE PMLPKETVLSLKRYNLGSSAMKRNVPGHVLQRPSYLTIRIQVTLLCNSSAEAL
Research Area	Cell Biology
Source	E.coli
Target Names	BLID
Protein Names	Recommended name: BH3-like motif-containing cell death inducer Alternative name(s): Breast cancer cell protein 2
Expression Region	1-108aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	39.0kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.