




Recombinant Human Galactose-3-O-sulfotransferase 3 (GAL3ST3), partial

Product Code	CSB-EP822172HU
Abbreviation	Recombinant Human GAL3ST3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96A11
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	SWYPKLFPLSCPPLRNSPPRPKHMTVAFLKTHKTAGTTVQNILFRFAERHNLTV VALPHPSCEHQFCYPRNFSAHFVHPATRPPHVLASHLRFDRAELERLMPPSTV YVTILREPAAMFESLFSYYNQYCPAFRRVNPASLEAFLRAPEAYYRAGEHFAM FAHNTLAYDLGGDNERSPRDDAAYLAGLIRQVEEVFSLVMIAEYFDESLVLLRR LLAWDLDDVLYAKLNARAASSRLAAIPAALARAARTWNALDAGLYDHFNATFW RHVARAGRACVEREARELREARQRLLRRCFGDEPLLRPAAQIRTKQLQPWQP SRKVDIMGYDLPGGGAGPATEACLKLAMPEVQYSNYLLRKQKRRGGARARPE PVLNPPPRPIRVLPRGPQGP
Research Area	Others
Source	E.coli
Target Names	GAL3ST3
Expression Region	41-431aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	51.9 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Description

In the process of expressing the recombinant Human GAL3ST3 protein, e.coli cells are transfected with a DNA expression vector that carries the gene encoding the desired protein (41-431aa). The cells are then cultured to induce the expression of the intended protein. The protein features a N-terminal 10xHis tag and C-terminal Myc tag fusion. The recombinant Human GAL3ST3 protein is collected and purified from the cell lysate through affinity purification, with its purity surpassing 85%, as assessed by SDS-PAGE.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.