




Recombinant Human Short transient receptor potential channel 4-associated protein (TRPC4AP)

Product Code	CSB-EP840581HU
Abbreviation	Recombinant Human TRPC4AP protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TEL6
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	AAAPVAAGSGAGRGRRSAATVAAWGGWGGRRPRPGNILLQLRQGQLTGRGL VRAVQFTETFLTERDKQSKWSGIPQLLLKLHTTSHLHSDFVEQCNILKEISPLLS MEAMAFVTEERKLTQETTYPNTYIFDLFGGVDLLVEILMRPTISIRGQKLKISDE MSKDCLSILYNTCVCTEGVTKRLAEKNDFVIFLFTLMTSKKTFLQTATLIEDILGV KKEMIRLDEVPNLSSLVSNFDQQQLANFCRILAVTISEMDTGND DKHTLLAKNA QQKKSLSLGPSSAAEINQAALLSIPGFVERLCKLATRKVSESTGTASFLQELEE YTWLDNALVLDALMRVANEESEHNQASIVFPPPGASEENGLPHTSARTQLPQ SMKIMHEIMYKLEVLYVLCVLLMGRQRNQVHRMIAEFKLIPGLNNLFDKLIWRK HSASALVLHGHNNQNCDCSPDITLKIQLRLLQSFSDDHNNKYLLNNQELNELS AISLKANIPVEEAVLNTDRSLVCDGKRGLLTRLLQVMKKEPAESSFRFWQARA VESFLRGTTSYADQMFLKRGLLLEHILYCIDVSECKSRDVLQSYFDLLGELMKF NVDAFKRFNKYINTDAKFQVFLKQINSSLVDSNMLVRCVTLSDRFENQVDMK VAEVLSECRLLAYISQVPTQMSFLFRLINIIHVQTLTQENVSCNLNTSLVILMLARR KERLPLYLRLLRMEHSKKYPGFLNNFHNLLRFWQQHYLHKDKDSTCLENS SCISFSYWKETVSILLNPDRQSPSALVSYIEEPPYMDIDRDFTEE
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	TRPC4AP
Expression Region	2-797aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged



Mol. Weight	94.8 kDa
Protein Length	Full Length of Mature Protein

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-----------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
-------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------