






Recombinant Human Terminal nucleotidyltransferase 4B (TENT4B)

Product Code	CSB-EP843329HU
Abbreviation	Recombinant Human TENT4B protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NDF8
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MYRSGERLLGSHALPAEQRDFLPLETTNNNNHHQPGAWARRAGSSASSPP SASSSPHPSAAVPAADPADSASGSSNKRKRDNKASGGRAAGGGRADGGGV VYSGTPWKRRNYNQGVVGLHEEISDFYEYMSPRPEEEKMRMEVVNRIESVIK ELWPSADVQIFGSFKTGLYLPTSDIDL VVFGKWENLPLWTL EEALRKHKVADE DSVKVLDKATVPIIKLTDSFTEVKVDISFNVQNGVRAADLIKDFTKKYPVLPYLV VLKQFLLQRDLNEVFTGGIGSYSLFMAVSFLQLHPREDACIPNTNYGVLLIEFF ELYGRHFNYLKTGIRIKDGGSYVAKDEVQKNMLDGYRPSMLYIEDPLQPGNDV GRSSYGAMQVKQAFDYAYVVL SHAVSPIAKYYPNNETESILGRIIRVTDEVATY RDWISKQWGLKNRPEPSCNGPVSSSSATQSSSSD VSDATPCKTPKQLLCRP STGNRVGSQDVSLESSQAVGKMQSTQTTNTSNSTNKSQHGSARLFRSSSKG FQGTQTSHGSLMTNKQH QHGKSNNQYYHGKKRKHKRDAPLSDL CR
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	PAPD5
Protein Names	Recommended name: PAP-associated domain-containing protein 5 EC= 2.7.7.- Alternative name(s): Terminal uridylyltransferase 3 Short name= TUTase 3 Topoisomerase-related function protein 4-2 Short name= TRF4-2
Expression Region	1-572aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	67.3 kDa
Protein Length	Full Length
Image	



The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.