




Recombinant Human E3 ubiquitin-protein ligase pellino homolog 3 (PELI3)

Product Code	CSB-EP847618HU
Abbreviation	Recombinant Human PELI3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N2H9
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MVLEGNPEVGSPRTSDLQHRGNKGSCVLSSPGEDAQPGEEPIKYGELIVLGC CEEGGEETEAQRGEVTGPRAHSCYNGCLASGDKGRRRSRLALSRRSHANGV KPDVMHHISTPLVSKALSNRGQHSISYTLRSHSVIVEYTHDSDTDMFQIGRST ENMIDFVVTDTSPGGGAAEGPSAQSTISRYACRILCDRRPPYTARIYAAGFDAS SNIFLGERAAKWRTPDGLMDGLTTNGVLVMHPAGGFSEDSAPGVWREISVCG NVYTLRDSRSAQQRGKLVENESNVLQDGLIDLCGATLLWRTPAGLLRAPTLK QLEAQRQEANAARPQCPVGLSTLAFPSARGRTAPDKQQPWVYVRCGHVHG YHGWGCRRERGPPQERECPLCRLVGPIVPLWLWGQEAGLCCLDPGPPSHAFAPC GHVCSEKTARYWAQTPLPHGTHAFHAACPFCAWLTGEHGCVRILFQGPLD
Research Area	Immunology
Source	E.coli
Target Names	PELI3
Protein Names	Recommended name: E3 ubiquitin-protein ligase pellino homolog 3 Short name= Pellino-3 EC= 6.3.2.-
Expression Region	1-469aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	54.8 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.