




Recombinant Human Charged multivesicular body protein 7 (CHMP7)

Product Code	CSB-EP855478HUa0
Abbreviation	Recombinant Human CHMP7 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WUX9
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MWSPEREAEAPAGGDPAGLLPPEWEEDEERMSFLFSAFKRSREVNSTDWDS KMGFWAPLVLSHSRRQGVRRLRLDLQEAQFQRKGSVPLGLATVLQDLLRRGE LQRESDFMASVDSSWISWGVGVFLLKPLKWTLSNMLGDNKVPAAEEVLVAVEL LKEKAAEEVYRLYQNSPLSSHPVVALSELSTLCANSCPDERTFYLVLLQLQKEKR VTVLEQNGEKIVKFARGPRAKVSPVNDVDVGYYQLMQSEQLLSRKVESLSQE AERCKEEARRACRAGKKQLALRLSLKAKQRTEKRIEALHAKLDTVQGILDRIYAS QTDQMFVNAYQAGVGALKLSMKDVTVEKAESLVDQIQELCDTQDEVSTLAG GVTNGLDFDSEELEKELDILLQDTTKEPLDLPDNPNNRHFTNSVNPNNRISDAEL EAELEKLSLSEGGGLVPSSKSPKRQLEPTLKPL
Research Area	Signal Transduction
Source	E.coli
Target Names	CHMP7
Expression Region	1-453aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	56.9 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Description

The process for producing recombinant Human CHMP7 protein in e.coli cells follows a sequence of steps: creating the expression vector that codes for the Human CHMP7 protein (1-453aa), introducing it into e.coli cells, cultivating and inducing of positive cells for protein expression, lysing cells, and conducting expression analysis. The resultant recombinant Human CHMP7 protein is purified from the cell lysate using affinity purification, with its purity exceeding 85%, as determined by SDS-PAGE.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.