




# Recombinant Rat E3 ubiquitin-protein ligase TRIM63 (Trim63)

<b>Product Code</b>	CSB-EP856332RA
<b>Abbreviation</b>	Recombinant Rat Trim63 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q91Z63
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	MDYKSGLPDGNAMENLEKQLICPICLEMFTKPVVILPCQHNLCKRCANDIFQA ANPYWTNRGGSVMSGGRFRCPSRHEVIMDRHGVYGLQRNLLVENIIDYK QECSSRPLQKGSHPMCKEHEDEKINIYCLTCEVPTCSLCKVFGAHQACEVAPL QSIFQGQKTELSNCISMLVAGNDRVQTIISQLEDSCRVTKENS HQVKEELSHKF DALYAILDEKKSELLQRITQE QEEKLDFIEALILQYREQLEKSTKL VETAIQSLDE PGGATFLLSAKPLIKSIVEASKGCQLGKTEQGFENMDYFTLNLEHIAEALRAIDF GTDEEEEFTEEEEEEDQEEGVSTEGHQ
<b>Research Area</b>	Cardiovascular
<b>Source</b>	E.coli
<b>Target Names</b>	Trim63
<b>Expression Region</b>	1-351aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	47.2 kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Description

The recombinant Rat Trim63 protein is made by following sequential steps: inserting a DNA fragment encoding the Rat Trim63 protein (1-351aa) into a plasmid vector, transferring this vector into e.coli cells, selecting the positive cells, culturing positive cells, inducing the protein expression, and collecting the Trim63 protein. The protein features a N-terminal 10xHis tag and C-terminal Myc tag fusion. The recombinant Trim63 protein is subjected to affinity purification and SDS-PAGE analysis. Its purity is over 90%.

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.