




Recombinant Human E3 ubiquitin-protein ligase pellino homolog 1 (PELI1)

Product Code	CSB-EP856932HUa0
Abbreviation	Recombinant Human PELI1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96FA3
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MFSPDQENHPSKAPVKYGELIVLGYNGLPNGDRGRKSRFALFKRPKANGV KPSTVHIACTPQAAKAISNKDQHSISYTLSTRAQTVVVEYTHDSNTDMFQIGRST ESPIDFVVTDTVPGSQSNSDTQSVQSTISRFAICRIERNPPFTARIYAAGFDSS KNIFLGEKAAKWKTSDGQMDGLTTNGVLVHPRNGFTEDSKPGIWREISVCG NVFSLRETRSAQQRGKMVEIETNQLQDGLIDLCGATLLWRTAEGLSHTPTVK HLEALRQEINAARPQCPVGFNTLAFPSMKRKDVVDEKQPWWVYLNCGHVHGYH NWGNKEERDGDRECPMCRSVGPYVPLWLGCAGFYVDAGPPTHAFSPCG HVCSEKTTAYWSQIPLPHGTHTFHAACPFCAHQLAGEQGYIRLIFQGPLD
Research Area	Signal Transduction
Source	E.coli
Target Names	PELI1
Expression Region	1-418aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	50.4 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.