



Recombinant Mouse Sarcosine dehydrogenase, mitochondrial (Sardh), partial


Product Code	CSB-EP857487MO
Abbreviation	Recombinant Mouse Sardh protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
Uniprot No.	Q99LB7
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Mouse Sarcosine dehydrogenase, mitochondrial (Sardh), partial
Immunogen Species	Mus musculus (Mouse)
Sensitivity	Not Test
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	VVVIGGGSLGCQTLYHLAKLGVGGA VLLERERLTSGTTWHTAGLLWQLRPSD VEVELLAHTRQVVSRDLEEETGLHTGWIQNGGLFIASNQQR LNEYKRLMSLGK AYGIESHVLSPAETKSLYPLMNVDDLYGTLYVPQDGTMDPAGTCTTLTRAAVA RGAQVIENCAVTGIRVRTDDFGVRRVA AVETE H GSIQTPCVVNCAGVWASKV GRMAGVKVPLVAMHHAYVVTERIEGIQNMPNVRDHDASVYLRLQGDALSVGG YEANPIFWEEVSDKFAFGLFDLDWDVFTQHIEGAINRVPVLEKTGIKSTVCGPE SFTPDHKPLMGEAPELRGFFLGCGFNSAGMMLGGGCGQELAHWIVHGRPEK DMYSYDIRRFHHS LTDHTRWIRERSHESYAKNYSVVFPHDEPLAGRNMRRDP LHEELLGQGCVFQERQGWERPWFNPQETAQVLDYDYYGAYGNQAHKDYT YSRLLGDEYTFDFPPHHHMIQKECLACRGAAAVFNMSYFGKFYLLGVDARKA ADWLFSADVNRPPGSTVYTCMLNQRGGTESDLTVSRLAPGTQASPLVPAFEG DCYYLAVGGGAVAQHNWSHINTVLQDQEFRCQLMDSSEDLGMLSIQGPASRDI LQDVLDADLSNEAFPSTHQLVRAAGHLVRAIRLSFVGELGWELHVPRASCLP VYRAVMAAGARHGLVNAGYRAIDSLSIEKGYRHWHADLRPDDSPLEAGLAF
Research Area	Metabolism
Source	E.coli
Target Names	Sardh
Expression Region	70-798aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	C-terminal 6xHis-tagged



Mol. Weight 87.7 kDa

Protein Length Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.