




Recombinant Human Kelch domain-containing protein 3 (KLHDC3)

Product Code	CSB-EP858414HUa0
Abbreviation	Recombinant Human KLHDC3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BQ90
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Kelch domain-containing protein 3(KLHDC3)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MLRWTVHLEGGPRRVNHA AVAVGHRVYSFGGYCSGEDYETLRQIDVHIFNAV SLRWTKLPPVKS AIRGQAPVVPYMRYGHSTVLIDDTVLLWGGRNDTEGACNV LYAFDVNTHKWFTPRVSGTVP GARDGHSACVLGKIMYIFGGYEQQADCFSNDI HKLDTSTMTWTLIC TKGSPARWRDFHSATMLGSHMYVFGGRADRF GPFHNSN NEIYCNIRIRVFDTRTEAWLDCPPTPVLPEGRRSHSAFGYNGELYIFGGYNARL NRHFHDLWKFNVPVSFTWKKIEPKGKGPCPRRRQCCCVGDKIVLFGGTSPSP EEGLGDEFDLIDHSD LHILDFSPSLKTLCKLAVIQYNLDQSCLPHDIRWELNAMT TNSNISRP IVSSHG
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	KLHDC3
Expression Region	1-382aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	49.1
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.