




Recombinant Human Hydroperoxide isomerase ALOXE3 (ALOXE3)

Product Code	CSB-EP863156HU
Abbreviation	Recombinant Human ALOXE3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BYJ1
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Hydroperoxide isomerase ALOXE3(ALOXE3)
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MAVYRLCVTTGPYLRAGTLDNISVTLVGTCGESPKQRLDRMGRDFAPGSVQK YKVRCTAELGELLLL RVHKERYAFFRKDSWYCSRICVTEPDGSVSHFPCYQWI EGYCTVELRPGTARTICQDSLPLLLDHRTRELRARQECYRWKIYAPGFPCMVD VNSFQEMESDKKFALTKTTTCVDQGDSSGNRYLPGFPMKIDIPSLMYMEPNV RYSATKTISLLFNAIPASLGMKLRGLLDRKGSWKKLDDMQNIFWCHKTFTTKYV TEHWCEDHFFGYQYLNQVNPVMLHCISLPSKLPVTNDMVAPLLGQDTCLQT ELERGNIFLADYWILAEAPTHCLNGRQQYVAAPLCLLWLSPQGALVPLAIQLSQ TPGPDSPIFLPTDSEWDWLLAKTWVRNSEFLVHENNTHFLCTHLLCEAFAMAT LRQLPLCHPIYKLLPHTRYTLQVNTIARATLLNPEGLVDQVTSIGRQGLIYLMS TGLAHFTYTNFCLPDSLRLRARGVLAIPNYHYRDDGLKIWAAIESFVSEIVGYYYPS DASVQQDSELQAWTGEIFAQAFGLGRESSGFPSRLCTPGEMVKFLTAIIFNCSA QHA AVNSGQHDFGAWMPNAPSSMRQPPPQTKGTTTLKTYLDTLPEVNI SCNN LLFWLV SQEPKDQRPLGTYPDEHFTEEAPRRSIAAFQSRLAQISRDIQERNQ GLALPYTYLDPPLIENSVSI
Research Area	Cell Biology
Source	E.coli
Target Names	ALOXE3
Expression Region	1-711aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	87.4 kDa


Protein Length

Full Length

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.