




Recombinant Pisum sativum Zinc finger CCCH domain-containing protein ZFN-like

Product Code	CSB-EP866006EWE
Abbreviation	Recombinant Pisum sativum Zinc finger CCCH domain-containing protein ZFN-like protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SWF9
Product Type	Recombinant Protein
Immunogen Species	Pisum sativum (Garden pea)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MEFDTIPQEAMWQMNLRSSETMESGPYPEHPGEPDCSYIIRTGLCRFGATCR FNHPPNRKLAIATARMKGEFPERLGQPECQYYLKTGTCKFGATCRFHHPKDK AGVAGRVALNILGYPLRPNESERAYYLRTGQCKFGNTCKFHHPQPSNMVLSM RGSTVYPTVQSPTTPGQQSYAAGITNWSSSSYVSPRWQGPSSYAPLILPQG VVSVPGWSTYGGQMGESEPQQTMRNDQTYGTSHQGELENAGLQGAYSQYR SGSVPVGFYALQRDNIFPERPDQPECQFYMKTDCKFGAVCRFHHPRERQIP APDCVLSPIGLPLRPGEPLCVFYSRYGICKFGPSCKFDHPMGIFTYNVASPLAD TPGRRLLGSSSGTAALSLSSEGLVESGTAKPRRLSLSETRPIPPGDDNIDDEG
Research Area	Others
Source	E.coli
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein ZFN-like
Expression Region	1-417aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	53.0 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.