



Recombinant Human CCR4-NOT transcription complex subunit 8 (CNOT8)

Product Code	CSB-EP871398HU
Relevance	Has 3'-5' poly(A) exoribonuclease activity for synthetic poly(A) RNA substrate. Its function seems to be partially redundant with that of CNOT7. Catalytic component of the CCR4-NOT complex which is linked to various cellular processes including bulk mRNA degradation, miRNA-mediated repression, translational repression during translational initiation and general transcription regulation. During miRNA-mediated repression the complex seems also to act as translational repressor during translational initiation. Additional complex functions may be a consequence of its influence on mRNA expression. Associates with members of the BTG family such as TOB1 and BTG2 and is required for their anti-proliferative activity.
Abbreviation	Recombinant Human CNOT8 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UFF9
Alias	CAF1-like protein
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MPAALVENSQVICEVWASNLEEEEMRKIREIVLSYSYIAMDTEFPGVVVRPIGEF RSSIDYQYQLLR CNVDLLKIIQLGLTFTNEKGEYPSGINTWQFNFKFNLTEDMY SQDSIDLLANSGLQFQKHEEEGIDTLHFAELLM TSGVVLCDNVKWLSFHSGYD FGYMKVLLTDSRLPEEEHEFFHILNLFPSIYDV KYLMKSCKNLKGGLQEVA DQ LDLQRIGRQH QAGSDSLTGM AFFRMKELFFEDS IDDAKYCGRLYGLGTGVA QKQNE DVDSAQEKMSILAIINMQQ
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	CNOT8
Protein Names	Recommended name: CCR4-NOT transcription complex subunit 8 Alternative name(s): CAF1-like protein Short name= CALIFp CAF2 CCR4-associated factor 8
Expression Region	1-292aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.




Tag Info N-terminal GST-tagged

Mol. Weight 60.5kDa

Protein Length Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.