




Recombinant Human Protein FAM49B (FAM49B)

Product Code	CSB-EP873658HU
Abbreviation	Recombinant Human FAM49B protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NUQ9
Alias	L1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MGNLLKVLCTDLEQGPNFFLDFENAQPTSESEKEIYNQVNVVLKDAEGILEDLQ SYRGAGHEIREAIQHHPADEKLQEKAAGVAVPLVGLKLFYEFSSQRLEAALRGL LGALTSTPYSPTQHLEREQALAKQFAEILHFTLRFDCLKMTNPAIQNDFSYYRR TLSRMRINNVPAGEGENEVNNELANRMSLFYAEATPMLKTLSDATTKFVSENKN LPIENTTDCLSTMASVCRVMLETPEYRSRFTNEETVSFCLRVMVGVIIYDHHVH PVGAFAKTSKIDMKGCIVLKDQPPNSVEGLLNALRYTTKHLNDETTSKQIKSM LQ
Research Area	Cell Biology
Source	E.coli
Target Names	FAM49B
Expression Region	1-324aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	52.7kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.