




Recombinant Human p53-regulated apoptosis-inducing protein 1 (TP53AIP1)

Product Code	CSB-EP884624HU
Relevance	May play an important role in mediating p53/TP53-dependent apoptosis.
Abbreviation	Recombinant Human TP53AIP1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HCN2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MGSSSEVSFRSAQASCSGARRQGLGRGDQNLSVMPPNGRAQTHTPGWVSP CSEN RDGLLPATAPGRLCSHRGADIPSFQTHQDPVTASGSSELHADCPQFRA LDRAGN
Research Area	Cell Biology
Source	E.coli
Target Names	TP53AIP1
Protein Names	Recommended name: p53-regulated apoptosis-inducing protein 1 Short name= p53AIP1
Expression Region	1-108aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	38.3kDa
Protein Length	Full Length of BC069399
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.