




# Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At4g32450, mitochondrial (PCMP-H63)

<b>Product Code</b>	CSB-EP890450DOA
<b>Abbreviation</b>	Recombinant Mouse-ear cress PCMP-H63 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9SUU7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	QNWQSSDGCSSYGTTGNGVPEENNTGGNHFQQDHSGHSSLDELDSICREG KVKKAVEIIKSWRNEGYYVDLPRLFWIAQLCGDAQALQEAKVVHEFITSSVGIS DISAYNSIIMYSGCGSVEDALTVFNSMPERNLETWCGVIRCFKNGQGDAI DTFSRFKQEGNKPDGEMFKEIFFACGVLGDMNEGLLHFESMYKEYGIIPCMEH YVSLVKMLAEPGYLDEALRFVESMEPNVDLWETLMNLSRVHGDILGDRCD MVEQLDASRLNKEKAGLVPVKSSDLVKEKLQRMAGPNYGIRYMAAGDISR PENRELYMALKSLKEHMIEIGYVPLSKLALHDVDQESKDENLFNHNERFAFIST FLDTPARSLIRVMKNLRVCADCHNALKLMSKIVGRELISRDAKRFHHMKDGV SCREYW
<b>Research Area</b>	others
<b>Source</b>	E.coli
<b>Target Names</b>	PCMP-H63
<b>Protein Names</b>	Recommended name: Pentatricopeptide repeat-containing protein At4g32450, mitochondrial
<b>Expression Region</b>	111-537aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	53.3kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.