




Recombinant Human ATP-dependent RNA helicase DDX19B (DDX19B)

Product Code	CSB-EP892352HUb1
Abbreviation	Recombinant Human DDX19B protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UMR2
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human ATP-dependent RNA helicase DDX19B(DDX19B)
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	ATDSWALAVDEQEAAEAESLSNLHLKEEKIKPDTNGAVVKTNANAECTDEEEKE DRAAQSLNKLIRSNLVDNTNQVEVLQRDPNSPLYSVKSFEELRLKPQLLQGV YAMGFNRPSKIQENALPLMLAEPQNLIAQSQSGTGKTA AFLVLA MLSQVEPAN KYPQCLCLSPTYELALQTGKVIEQMGKFYPELKLAYAVRGNKLERGQKISEQIV IGTPGTVLDWCSKLKFIDPKKIKVFVLDEADVMIA TQGHQDQSIRIQRMLPRNC QMLLFSATFEDSVWKFAQKVVPDPNVIKLKREEETLDTIKQYYVLCSSRDEKF QALCNLYGAITIAQAMIFCHTRKTASWLAAELSKEGHQVALLSGEMMVEQRAA VIERFREGKEKVLVTTNVCARGIDVEQVSVVINFDLPVDKDG NPDNETYLHRIG RTGRFGKRGLAVNMVDSKHS MNILNRIQE HFNKKIERLDTDDLDEIEKIAN
Research Area	Others
Source	E.coli
Target Names	DDX19B
Expression Region	2-479aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	61.2 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.