




Recombinant Human Zinc finger and BTB domain-containing protein 32 (ZBTB32)

Product Code	CSB-EP897538HU
Abbreviation	Recombinant Human ZBTB32 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y2Y4
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MSLPPIRLPSPYGS DRLVQLAARLRPALCDTLITVGSQEFP AHS LVL AGVSQQ L GRRGQWALGEGISPSTFAQLLN FVYGESVELQPGELRPLQEAA RALGVQSLE EACWRARGDRAKKPD PGLKKHQEEPEKPSRNPERELGDPGEKQKPEQVSRT GGREQEMLHKHSPPRGRPEMAGATQEAQQEQTRSKEKRLQAPVGQRGADG KHGVLTLWLRNPGGSEESLRKLPGLPPAGSLQTSVTPRPSWAEAPWL VGG QPALWSILLMPPRYGIPFYHSTPTTGAWQEVWREQRIPLSLNAPKGLWSQNQ LASSSPTPGSLPQGP AQLSPGEMEESDQGHTGALATCAGHEDKAGCPPRPH PPPAPPARSRPYACSVCGKRFS LKHQMETHYRVHTGEKPFSCSLCPQRSRDF SAMTKHLRTHGAAPYRCSLCGAGCPSLASMQA HMRGHSPSQLPPGWTIRST FLYSSSRPSRPSTSPCCPSSSTT
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	ZBTB32
Protein Names	Recommended name: Zinc finger and BTB domain-containing protein 32 Alternative name(s): FANCC-interacting protein Fanconi anemia zinc finger protein Testis zinc finger protein Zinc finger protein 538
Expression Region	1-487aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	60.4 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.