




Recombinant Mouse Apolipoprotein E (ApoE)

Product Code	CSB-MP001936MO
Abbreviation	Recombinant Mouse ApoE protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08226
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	EGEPEVTDQLEWQSNQPWEQALNRFWDYLRWVQTLSDQVQEELQSSQVTQ ELTALMEDTMTEVKAYKKELEEQLGPVAAETRARLGKEVQAAQARLGADMED LRNRLGQYRNEVHTMLGQSTEEIRARLSTHLRKMRLMRDAEDLQKRLAVY KAGAREGAERGVSAIRERLGPLVEQGRQRTANLGAGAAQPLRDRAQAFGDRI RGRLEEVGNQARDRLEEVRHEMEEVRSKMEEQTQQIRLQAEIFQARLKGWFE PIVEDMHRQWANLMEKIQASVATNPITPVAQENQ
Research Area	Cardiovascular
Source	Mammalian cell
Target Names	ApoE
Protein Names	Recommended name: Apolipoprotein E Short name= Apo-E
Expression Region	19-311aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	39
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.