




Recombinant Human Desmoglein-3 (DSG3), partial,Biotinylated

Product Code	CSB-MP007205HUj9-B
Abbreviation	Recombinant Human DSG3 protein, partial, Biotinylated
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32926
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Desmoglein-3(DSG3),partial,Biotinylated
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	EWVKFAKPCREGEDNSKRNP IAKITSYQATQKITYRISGV GIDQPPFGIFVVDKNTG DINITAIVDREETPSFLIT CRALNAQGLDVEKPLIL TVKILDINDNPPVFSQQ IFMGEIEENSASNSLVM ILNATDADEPNHLNSKIA FKIVSQEPAGTPMFLLS RNTGEVRTLTNSLDREQ ASSYRLVVGADKDGEGL STQCECNKVKDVNDNF PMFRDSQYSARIEENIL SSELLRFQVTDLDEEYTD NWLAVYFFTSGNEGNWF EIQTDPRTNELKVVKA LDYEQLSVKLSIAVKNK AEFHQSVISRYRVQSTP VTIQVINVREGIAFRP ASKTFTVQKGISSKKLV DYLGTQAIDEDTNKAAS NVKYVMGRNDGGYLM IDSKTAEIKFVKNMNRD STFIVNKTITAEVLAID EYTGKTSTGTVYVRVP DFNDNCPTAVLEKDAVC SSSPSVVVSARTLNNRY TGPYTFALEDQPVKLPA VWSITTLNATSALLRAQ EQIPPGVYHISLVLTDS QNNRCEMPLSLTLEVCQ CDNRGICGTSYPTTSPG TRYGRPHSGR
Research Area	Cell Adhesion
Source	Mammalian cell
Target Names	DSG3
Expression Region	50-615aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-Avi-tagged
Mol. Weight	67.5 kDa
Protein Length	Partial



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.