




Recombinant Mouse GDNF family receptor alpha-1 (Gfra1)

Product Code	CSB-MP009379MO
Abbreviation	Recombinant Mouse Gfra1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P97785
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	DRLDCVKASDQCLKEQSCSTKYRRTLRCQCVAGKETNFSLTSGLEAKDECRSAM EALKQKSLYNCRCKRGMKKEKNCLRIYWSMYQSLQGNDLLEDSPYEPVNSRL SDIFRAVPFISDVFQQVEHISKGNCLDAKACNLDDTCKKYRSAYITPCTTSM SNEVCNRRKCHKALRQFFDKVPAKHSYGMLFCSCRDVACTERRRQTIVPVCS YEERERPNCNLQDSCKTNYICRSRLADFFTNCQPESRSVSNCLKENYADCLL AYSGLIGHTVMTPNYIDSSSLSVAPWCDCSNSGNDLEDCLKFLNFFKDNTCLKN AIQAFGNGSDVTMWQPAPPVQTTTATTTTAFRIKNKPLGPAGSENEIPTHVLPP CANLQAQKLKSNVSGSTHLCLSDNDYGKDGLAGASS
Research Area	Neuroscience
Source	Mammalian cell
Target Names	Gfra1
Protein Names	Recommended name: GDNF family receptor alpha-1 Short name= GDNF receptor alpha-1 Short name= GDNFR-alpha-1 Short name= GFR-alpha-1 Alternative name(s): TGF-beta-related neurotrophic factor receptor 1
Expression Region	25-430aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	47.6
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.