





Recombinant Mouse Inositol 1,4,5-trisphosphate receptor-interacting protein-like 1 (Itpril1), partial (Active)

Product Code	CSB-MP011916MO1
Abbreviation	Recombinant Mouse Itpril1 protein, partial (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2ASA8
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, 6% Trehalose, pH 7.4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Biological Activity	Measured by its binding ability in a functional ELISA. Immobilized Mouse Itpril1 at 2 µg/mL can bind Anti-ITPRIPL1 recombinant antibody (CSB-RA756966MA1HU). The EC50 is 0.7534-0.9877 ng/mL.
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	DRMDLDTLARSRLQLEKRMSEEMRQLEMEFEERSRAAEQKQKVENFWRGDTS SDQLVLGKKDMGWPFQAGGQDGG
Source	Mammalian cell
Target Names	Itpril1
Expression Region	23-96aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	10.0 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Activity

Measured by its binding ability in a functional ELISA. Immobilized Mouse Itpril1 at 2 µg/ml can bind Anti-ITPRIPL1 recombinant antibody (CSB-RA756966MA1HU). The EC₅₀ is 0.7534-0.9877 ng/mL.

Endotoxin

Less than 1.0 EU/ug as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.