



# Recombinant Human Vascular endothelial growth factor receptor 2 (KDR), partial

<b>Product Code</b>	CSB-MP012145HU
<b>Abbreviation</b>	Recombinant Human KDR protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P35968
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
<b>Product Type</b>	Recombinant Human Vascular endothelial growth factor receptor 2 (KDR), partial
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sensitivity</b>	Not Test
<b>Purity</b>	Greater than 95% as determined by SDS-PAGE. Greater than 90% as determined by SEC-HPLC.
<b>Sequence</b>	ASVGLPSVSLDLPRLSIQKDILTIKANTTLQITCRGQRDLWLWPNNQSGSEQR VEVTECSDGLFCKTLTIPKVIGNDTGAYKCFYRETDLASVIYVYVQDYRSPFIAS VSDQHGVVYITENKNKTVVIPCLGSISNLNVSLCARYPEKRFVDPGNRISWDSK KGFTIPSYMISYAGMVFCEAKINDESYQSIMYIVVVVGYRIYDVVLSPSHGIELS VGEKLVLNCTARTELVGIDFNWEYPSSKHQHKLVNRDLKTQSGSEMKKFL STLTIDGVTRSDQGLYTCAASSGLMTKKNSTFVRVHEKPFVAFGSGMESLVEA TVGERVRIPAKYLGYPPEIKWYKNGIPLESNHTIKAGHVLTIMEVSERDTGNY TVILTNPISKEKQSHVVSLVYVPPQIGEKSLISPVDSYQYGTQTTLTCTVYAIPP PHIIHWYWQLEEECANEPSQAVSVTNPYPCEEWRSVEDFQGGNKIEVNKNQ FALIEGKNKTVSTLVIQAANVSALYKCEAVNKKVGRGERVISFHVTRGPEITLQPD MQPTEQESVSLWCTADRSTFENLTWYKLGQPQLPIHVGEPLTPVCKNLDTLW KLNATMFSNSTNDILIMELKNASLQDQGDYVCLAQDRKTKKRHCVVRQLTVLE RVAPTITGNLENQTTSIGESIEVSCTASGNPPPQIMWFKDNETLVEDSGIVLKD GNRNLTIRRVKEDGLYTCQACSVLGCAKVEAFFIIEGAQEKTNLE
<b>Research Area</b>	Cancer
<b>Source</b>	Mammalian cell
<b>Target Names</b>	KDR
<b>Expression Region</b>	20-764aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.




**Tag Info** C-terminal hFc1-tagged


**Mol. Weight** 112.1 kDa

**Protein Length** Partial

### Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



The purity of KDR was greater than 90% as determined by SEC-HPLC

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.