




Recombinant Human Putative protein SFTA3 (SFTA3)

Product Code	CSB-MP021167HU
Abbreviation	Recombinant Human SFTA3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C7M3
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MRAGFSDFQLIRDQVLFLQDQAQRLTEWLQLSGFENPVSESTTLCLREREKRI PTCVAVCVSPSGTVHTALLHPTTLSQSRSSSEAKMLIIHTA
Research Area	Cell Biology
Source	Mammalian cell
Target Names	SFTA3
Expression Region	1-94aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-Myc-tagged
Mol. Weight	14.5kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.