




Recombinant Macaca fascicularis Heat-stable enterotoxin receptor (GUCY2C), partial

Product Code	CSB-MP3779MOV
Abbreviation	Recombinant Cynomolgus monkey GUCY2C protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	XP_005570270
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Macaca fascicularis Heat-stable enterotoxin receptor (GUCY2C), partial
Immunogen Species	Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey)?????macaque??????
Sensitivity	Not Test
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	SQVSNCHNGSYEISVLMMDNSAF AEPLENVEDAVNEGLEIVRGR LQNAGLN VTVNASFMYS DGLIHNSGDCRSSTCEGLDLLRKISNAKRMGCVLMGPSCTYST FQMYLDTELSYPMISAGSFGLSCDYKETLTRLMS PARKLTYFLVNFWKTN DLP FKTYSWSTSYVYKNGTESEDCFWYLNAL EASVSYSFHEL SFKLVL RQDKEFQ DILMDHNRKSNVIVMCGDPQFLYKLKGDR AAVEDIIVLVDLFNDQYFEDN VTA PDYMKNVLVLTQSPGN SLLNSSFSRNLSPTKRDFALAYLNGILLFGHMLKTFLE NGENITTPKFAHA FRNLTFEGYDGPVTLDDWGDVDSTMVLLYTSVDTKKYKVL LTYDTHVNQTNPVDMSP TFTWKNSKLPNDITDRGPQ
Research Area	Others
Source	Mammalian cell
Target Names	GUCY2C
Expression Region	24-430aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	48.3 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.