





# Recombinant Macaca fascicularis Leukemia inhibitory factor (LIF)

|                          |                                                                                                                                                                                                                                                                                 |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-MP5062MOV                                                                                                                                                                                                                                                                   |
| <b>Abbreviation</b>      | Recombinant Cynomolgus monkey LIF protein                                                                                                                                                                                                                                       |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| <b>Uniprot No.</b>       | A0A2K5UCA8                                                                                                                                                                                                                                                                      |
| <b>Form</b>              | Liquid or Lyophilized powder                                                                                                                                                                                                                                                    |
| <b>Storage Buffer</b>    | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.                                                       |
| <b>Product Type</b>      | Recombinant Macaca fascicularis Leukemia inhibitory factor (LIF)                                                                                                                                                                                                                |
| <b>Immunogen Species</b> | Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey)                                                                                                                                                                                                                   |
| <b>Sensitivity</b>       | Not Test                                                                                                                                                                                                                                                                        |
| <b>Purity</b>            | Greater than 95% as determined by SDS-PAGE.<br>Greater than 95% as determined by SEC-HPLC.                                                                                                                                                                                      |
| <b>Sequence</b>          | SPLPITPVNATCAIRHPCHNNLMNQIRSQLAQLNGSANALFILYYTAQGEPFPN<br>NLDKLCGPNVTDFFPPFHANGTEKAKLVELYRIVVYLGTS LGNITRDQKILNPSAL<br>SLH SKLNATADILRGLLSNVLCRLCSKYHVGHV DVTYGPDTSGKDVFQKKKLG<br>CQLLGKYKETIAVLAQAF                                                                            |
| <b>Research Area</b>     | Others                                                                                                                                                                                                                                                                          |
| <b>Source</b>            | Mammalian cell                                                                                                                                                                                                                                                                  |
| <b>Target Names</b>      | LIF                                                                                                                                                                                                                                                                             |
| <b>Expression Region</b> | 23-202aa                                                                                                                                                                                                                                                                        |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                             |
| <b>Tag Info</b>          | C-terminal 10xHis-tagged                                                                                                                                                                                                                                                        |
| <b>Mol. Weight</b>       | 23.4 kDa                                                                                                                                                                                                                                                                        |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                   |
| <b>Image</b>             |                                                                                                                                                                                                                                                                                 |



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



The purity of LIF was greater than 95% as determined by SEC-HPLC

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.