




Recombinant Human tRNA methyltransferase 10 homolog C (TRMT10C)

Product Code	CSB-MP748718HU
Abbreviation	Recombinant Human TRMT10C protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7L0Y3
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MSSKIPAVTYPKNESTPPSEELELDKWKTTMKSSVQEECVSTISSSKDEDPLA ATREFIEMWRLLGREVPHEHITEEELKTLMECVSNTAKKKYLKYLTYKEKVKKAR QIKKEMKAAAREEAKNIKLLETTEEDKQKNFLFLRLWDRNMDIAMGWKGAQA MQFGQPLVFDMAYENYMKRKELQNTVSQELLESEGWNRRNVDPFHIYFCNLKI DGALHRELVKRYQEKWDKLLLTSTEKSHVDLFPKDSIIYLTADSPNVMTTFRHD KVVYIGSFVDKSMQPGTSLAKAKRLNLATECLPLDKYLQWEIGNKNLTLDQMIR ILLCLKNNGNWQEALQFVPRKHTGFLEISQHSQEFINRLKKAKT
Research Area	Others
Source	Mammalian cell
Target Names	TRMT10C
Protein Names	Recommended name: Mitochondrial ribonuclease P protein 1 Short name= Mitochondrial RNase P protein 1Alternative name(s): HBV pre-S2 trans-regulated protein 2 RNA (guanine-9-)-methyltransferase domain-containing protein 1 EC= 2.1.1.-
Expression Region	40-403aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	47.6 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.