




Recombinant Rat Amphoterin-induced protein 3 (Amigo3)

Product Code	CSB-MP802260RA
Relevance	May mediate heterophilic cell-cell interaction. May contribute to signal transduction through its intracellular domain.
Abbreviation	Recombinant Rat Amigo3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q80ZD5
Product Type	Recombinant Proteins
Immunogen Species	Rattus norvegicus (Rat)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	TSDLEGVLPPDPHNCPNKCVCAADVLSCAGRGLQDLPAALPATAAELDLSHN ALKRLHPGWLAPLSRLRALYLGYNKLDVLGRGVFTNASGLRILDLSNLLRRLR TYDLDGLEEELEKLLLFNNRLMHLDLDAFQGLSMLSHLYLSCNELSSFSFNHLH GLGLTRLRLTDLSSNWLGHVSVPELAALPTFLKNRLYLHNNPLPCDCSLYHLLR RWHQRGLSALHDFEREYTCIAFKVAESRVRFFEHSRVFKNCVAAAPGLELP EEELHTHVGGSLRLFCNTTVPAARVAWVSPKNELLVAPGSQDGSIAVLADGSL AIGRVQEQHAGLFVCLASGPRLHHNQTL EYNVSVHKPRPEPEAFNT
Research Area	Cell Biology
Source	Mammalian cell
Target Names	Amigo3
Expression Region	20-383aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-Myc-tagged
Mol. Weight	44.4kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.